

Empty Mixing Containers Accessories



STERILE EO

DMP Offers a Comprehensive Line of Flexible Empty Mixing Containers with:

- Sterile
- Large Bore (high flow)
- Flexible
- Latex-Free & DEHP-Free
- DPT6820 is a standard ISO Luer Fitment


Code	Units/ CS
DPT6820	50
DPT1130	30
DPT1133	50
DPT6720	50



DPT1133 - 36" 3-Lead Male Connector Set with Vented Spikes



DPT6820 – Port Male-to-Female Connector



DPT1130 -12" Extended Port Connector with Spike



DPT6720 - Port Female-to-Female Connector

Order: +1 844-323-9397
Product Info: +1 847-909-4298
+1 312-544-9691

www.ecoflx.com
sales@ecoflx.com