

# EcoFLX™

EMPTY MIXING CONTAINERS

Manufactured by Douglas Medical Products, Inc., IL, USA

## DML – EVA Light Sensitive Container:

- 3-Port Configuration
- Sterile
- Tamper-Evident Medication Site
- EVA Material (PVC-Free)
- Clear & Yellow & Flexible
- Latex-Free & DEHP-Free
- Administration Port
- Meet USP <671> Standard

Series: DML

Size mL	EVA Code	Units/ CS
100	DML0108*	50
250	DML0258*	50
500	DML0508*	50
1000	DML1008*	50
2000	DML2008*	30
3000	DML3008*	30

Empty  
Mixing  
Containers

STERILE R

Rx only

- *Note All sizes are Single Pack.*
- *DML0108 & DML3008 are custom order.*

[www.ecoflx.com](http://www.ecoflx.com)

[sales@ecoflx.com](mailto:sales@ecoflx.com)

Order: +1 844-323-9397  
Product Info: +1 847-909-4298

